



I – Identification of the Substance and of the Company	
<p>SUPPLIER: RMO, Inc. 650 W. Colfax Ave. Denver, CO 80204 303-592-8200</p> <p>Emergency Information Chemtrec: 800-424-9300 Chemtrec International: 202-483-7616</p> <hr/> <p>Product Grade / Name: Wax Base Product Code: N/A Chemical Family : MICROCRYSTALLINE PARAPHIN Chemical Name: FORMULATED PRODUCT Formula: N/A</p>	<p>Trade Name and Synonyms – Patient Wax</p> <p>Description: Wax Base</p>

II – Composition / Information on Ingredients			
Ingredient Name:			
Hazardous Ingredients: None			
/CAS Number	Exposure Limits	Concentration	
N/A	N/A	N/A	
Hazardous Materials Identification System (HMIS) (For Material as Packaged):			
HMIS RATINGS:	HEALTH	FLAMMABILITY	REACTIVITY
	0	0	0
0=MIMIMAL 1=SLIGHT 2=MODERATE 3=SERIOUS 4=SERIOUS			

III – Hazards Identification
Hazardous ingredients: This product has non-hazardous ingredient as defined under the criteria of the Federal OSHA Hazard Communication Standard 29 CFR 1911 1200
Medical Conditions
Aggravated by Exposure: None
Exposure Limits: None
Potential Health Effects: None
Route(s) of Entry: Ingestion
Human Effects and Symptoms of Overexposure: None
Acute Inhalation: N/A
Chronic Inhalation: None
Acute Eye Contact: N/A
Chronic Eye Contact: None

Carcinogenicity: None

IV – First Aid Measures

First aid for eyes: In case of contact, flush eyes well with water for 15 minutes. Obtain medical attention if irritation occurs.

First aid for Skin: Remove excess material from the skin with a waterless skin cleaner. Flush skin with plenty of water and soap.

First aid for inhalation: N/A

First aid for ingestion: Due to the physical state of this material, ingestion is unlikely to occur.

V – Fire Fighting Measures

FLASH POINT: Not combustible

FLAMMABLE LIMITS: N/A

Upper Explosive Limit (UEL): N/A

Over Explosive Limit (LEL): N/A

FIRE EXTINGUISHING MEDIA: Foam, Dry chemicals, CO₂, Water Fog, Sand/Earth

SPECIAL FIRE FIGHTING PROCEDURES: N/A

VI – Accidental Release Measures

Spill or leak procedures: Remove by mechanical means

VII – Handling and Storage

Storage temperature (min/Max: Ambient

Shelf Life: None

Special sensitive: Do not store or mix with strong oxidants

Handling/Storage Precautions: Store in sealed containers. Avoid excessive heat or direct sunlight.

VIII – Exposure Controls / Personal Protection

Eye Protection requirements: None

Skin Protection requirements: none required

Ventilation requirements: None

Respirator requirements: None

Additional protective measures: None

IX – Physical and Chemical Properties

Physical Form: Paste

Color	White	Orange	Yellow	Purple	Green	Red
Odor	Odorless	Orange	Lemon	Grape	Lime	Cherry

Ph: N/A

Boiling Point: N/A

Melting-Freezing point: 140/150F (Melting)

Solubility in water: Insoluble

Specific Gravity: 35/36G/CM³ (in Water=1)

Bulk Density: N/A

% Volatile by Volume: N/A

Vapor Pressure: N/A

X – Stability and Reactivity
Stability: Stable at ambient temperature and atmospheric pressure.
Hazardous Polymerization: Will not occur
Incompatibilities: None Known
Decomposition Products: By fire; carbon monoxide, carbon Dioxide, and various hydrocarbon fragments.

XI – Toxicological Information.
Toxicological testing has not been conducted with this material

XII – Ecological Information
None

XIII – Disposal Considerations
Waste Disposal (Insure Conformity with all Applicable Disposal Regulations): Incinerate or place in waste management facility.

XIV – Transportation Information
Technical Shipping Name: N/A
Freight Class Bulk: N/A
Freight Class Package: N/A
Product Label: N/A
Hazard Class or Division: Non-Hazardous
Hazardous Class Division Number: Not Hazardous by D.O.T. Regulations

XV – Regulatory Information
SARA Title III
Section 302 Extremely Hazardous Substances: None
Section 311/312 Hazardous Categories: None
Section 313 Toxic Chemicals: None
RERA Status: Non-Hazardous

XVI – Other Information
Note: While the information and recommendations set forth on this data sheet are believed to be accurate as received from our suppliers, RMO, Inc. makes no warranty with respect thereto and disclaims all liability from reliance thereon.